

[todo]+ DRIVE ORDER FORM

Date:

*Save this form and send it to sales@todo-works.com to submit your order.
**Orange, underlined text is [clickable](#) to show info and image.

Invoice Address		Shipping	
Company name		Ship to invoice address	Ship to alternate address
Street		Alternate Address	
		Street	
City	Zip		
Contact name	Country		
Phone		City	Zip
Email		Country	Remarks

Make your order

Item	Code No.	Description	Quantity
Todo Drive V2.1	TDV-V2.1	The world's lightest and smallest Power Add-on IoT device	
Todo Battery Pack	B010000	Lithium-ion Battery Pack, 144Wh, Solely for TodoDrive	



Click [this link](#) to see 'part list' to make extra orders.

Extra Order

Item	Code No.	Quantity

Your Wheelchair Information

Model			
Manufacture			

Image	Installation Criteria
	✓ Motor Mount Space must be min <u>distance 9cm</u> or above.
	✓ Todo Drive must not be mounted on <u>Carbon frame</u> wheelchairs. ✓ <u>Titanium frame</u> wheelchair may be damaged if Todo is installed.
	✓ Not ellipse-shaped frame but <u>circular-shaped frame</u>
	✓ <u>Camber</u> shall be no more than <u>3 degrees</u>.
	✓ Do not install Todo Drive on rear wheels below 20 inch. <u>Only installable for 20~26 inches rear wheels</u>

**Orange, underlined text is [clickable](#) to show its info and image.

[return to the order page](#)

Category	Code #	Item	MOQ
Main Components	M01M0000-N2	Motor Unit R	1
	M02M0000-N2	Motor Unit L	1
	C01C0000	Main Unit	1
	J01J0000	Joystick Unit	1
	B010000	Todo Battery-Pack, 144Wh	10
Motor	OP0004-1	Wheel 50mm AL	2
	M01M0022	Wheel 60mm PVC	2
	OP0001-1	Motor Body Pipe 65mm	2
	M01M0014-1	Motor Body Pipe 80mm	2
	OP0002-1	Motor Body Pipe 100mm	2
	OP0003-1	Motor Body Pipe 135mm	2
	M01M0016-1	Body Bracket(Normal)	2
	M01M0017-1	Upper Body Bracket 100mm	2
	OP0005-1	Upper Body Bracket 120mm	2
	OP0006-1	Upper Body Bracket(Diagonal)	2
	M01M0023-1	Lever 25	2
	M01M0024-1	Lever 40	2
	M01M0034	Motor Y-Cable	1
	M01M0032	Motor Andicap	20
	M01M0031	Motor Rubber Cover	2
	M01M0033	Adhesive Grip Strips	20
	OP0009	Wheel lock, standard, compl.	2
Joystick	ETC0019	J-Handle	10
	ETC0023	J-Gaiter Protector	10
	J01M0000-P	J-folding mount	1
	J01M0000-H	J-folding mount-2	1
Adapter	B010009	Todo Charger 16.8V 3.5Ah	1
Clamp	UC013555-2	W-Clamp S	2
	UC013565-2	W-Clamp M	2
	UC013585-2	W-Clamp L	2
	CI011734	W-Inner 17mm	10
	CI012234	W-Inner 22mm	10
	CI012534	W-Inner 25mm	10
	UC012065-2	T-Clamp S	2
	UC012085-2	T-Clamp L	2
	CI011719	T-Inner 17mm	10
	CI012019	T-Inner 20mm	10
	CI012219	T-Inner 22mm	10
	CI012519	T-Inner 25mm	10
	UC012066-1	Motus Clamp 20*70	1





Todo Drive | TDV-2.1

The world's lightest wheelchair Power Add-on IoT Device
The only 5kg kit attaching on your wheelchair gives you freedom of mobility.



+82-31-313-0351



www.todo-works.com



todo@todo-works.com



42- B-dong, Eungyejungang-ro 306beon-gil, 14923 Siheung-si, Republic of Korea



[todo]
WORKS



Todo Battery | B010000

Lithium-ion Battery pack, 14.4V X 10Ah = 144Wh



+82-31-313-0351



www.todo-works.com



todo@todo-works.com



42- B-dong, Eungyejungang-ro 306beon-gil, 14923 Siheung-si, Republic of Korea



[todo]
WORKS



Motor Unit | M01M000-N2(R), M02M000-N2(L)

138.5W DC Geared Motor



+82-31-313-0351



www.todo-works.com



todo@todo-works.com



42- B-dong, Eungyejungang-ro 306beon-gil, 14923 Siheung-si, Republic of Korea



[todo]
WORKS



1. J-Handle | ETC0019

2. J-Gaiter Protector | ETC0023

3. Joystick Unit | J01J0000





Main Unit | C01C0000



+82-31-313-0351



www.todo-works.com



todo@todo-works.com



42- B-dong, Eungyejungang-ro 306beon-gil, 14923 Siheung-si, Republic of Korea



[todo]
WORKS



1. Wheel 60mm PVC | M01M0022

2.Wheel 50mm AL | OP0-004-1





1. Motor Body Pipe 65mm | OP001-1

2. Motor Body Pipe 80mm | M01M0014-1

3. Motor Body Pipe 100mm | OP002-1

4. Motor Body Pipe 135mm | OP003-1





1. Body Bracket(Normal) | M01M0016-1

2. Upper Body Bracket(Diagonal) | OP0006-1

3. Upper Body Bracket 100mm | M01M0017-1

4. Motor Body Pipe 135mm | OP005-1



+82-31-313-0351



www.todo-works.com



todo@todo-works.com



42- B-dong, Eungyejungang-ro 306beon-gil, 14923 Siheung-si, Republic of Korea



**[todo]
WORKS**




- 1. Lever 25mm | M01M0023-1
- 2. Lever 40mm | MO1M0024-1
- 3. Motor Andicap | M01M0032

**Motor Y-Cable | M01M0034**



Motor Rubber Cover | M01M0031



Adhesive Grip Strips | M01C011



Wheel lock, standard, compl. | OP0009





J-Folding Mount | J01M000-V





J-Folding Mount 2 | J01M0000-H



+82-31-313-0351



www.todo-works.com



todo@todo-works.com



42- B-dong, Eungyejungang-ro 306beon-gil, 14923 Siheung-si, Republic of Korea





Todo Charger 16.8V 3.5Ah | B010009





1. W-Clamp S 35x55 | UC013555-2

2. W-Clamp M 35x65 | UC013565-2

3. W-Clamp L 35x85 | UC013585-2



+82-31-313-0351



www.todo-works.com




todo@todo-works.com



42- B-dong, Eungyejungang-ro 306beon-gil, 14923 Siheung-si, Republic of Korea



[todo]
WORKS



1. W-Inner 17mm | CI011734

2. W-Inner 22mm | CI012234

3. W-Inner 25mm | CI012534






1. T-Clamp S 20x65 | UC012065-2

2. T-Clamp L 20x85 | UC012085-2





1. T-Inner 22mm | CI012219

2. T-Inner 17mm | CI011719

3. T-Inner 20mm | CI012019

4. T-Inner 25mm | CI012519



+82-31-313-0351



www.todo-works.com



todo@todo-works.com



42- B-dong, Eungyejungang-ro 306beon-gil, 14923 Siheung-si, Republic of Korea



**[todo]
WORKS**



Motus Clamp 20x70 | UC012066-1

Applied to this clamp to mount Joystick unit in unique wheelchairs



+82-31-313-0351



www.todo-works.com



todo@todo-works.com



42- B-dong, Eungyejungang-ro 306beon-gil, 14923 Siheung-si, Republic of Korea



[todo]
WORKS